



lifeSign  
THE DXPRESS READER





## DXpress® Reader

*A multifunctional reader with user friendly features to deliver fast and accurate results in the laboratory or at the point-of-care*

- Easy to use; includes a barcode scanner for easy patient I.D. input
  - LOT specific test device information on data chip
  - Operator & QC lockout features
  - QC controls built into the reader software & test device
  - Stores & manages up to 1,000 patient and 60 QC results
- 
- On-board printer for hardcopy results
  - RS-232 port for LIS connectivity when downloading patient results
  - Instrument read eliminates visual subjectivity
  - Lightweight for portability, less than 3 lbs
  - Single or batch testing with on-board timer



### ORDERING INFORMATION

#### PRODUCT

#### ITEM #

DXpress® Reader

LSR-2000 (1 year warranty)

DXpress® Reader

LSR-2003 (3 year warranty)

## StatusFirst® CHF NT-proBNP

*For the rapid quantitative detection of NT-proBNP aiding in the diagnosis of congestive heart failure*



- Quantitative results in 15 minutes with DXpress® Reader
- Easy to use procedure; Add 3 drops of specimen into device
- EDTA plasma sample
- Measures a range of 20 to 5000 pg/mL
- Analytical sensitivity of 20 pg/mL

### ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® CHF NT-proBNP	20	60401
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® CHF NT-proBNP Controls	2 x 1 mL	LSCHF-300C

## StatusFirst® DOA

*One step immunoassay for the qualitative detection of drugs of abuse and their metabolites in urine*

- FDA cleared
- Simple one-step procedure
- Instrument read eliminates subjectivity
- Accurate results in 5 minutes
- Available in 5 & 10 drug panels
- Results stored and printed for documentation
- Walk away procedure accommodates busy labs and workflow interruptions



### ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® DOA 5 Panel THC/COC/OPI/MET/PCP	25	30525
StatusFirst® DOA 10 Panel MET/OPI/COC/THC/PCP/BZO/BAR/MTD/TCA/AMP	25	30025
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® DOA 10 Controls	2 x 5 mL	21000

## StatusFirst® hCG

*For the detection of hCG in serum or urine*



- Urine or Serum
- Superior sensitivity
- Simple procedure
- Instrument read eliminates subjectivity
- hCG results in 3-5 minutes

### ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® hCG Urine/Serum	35	35335
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® hCG Urine Controls	2 x 5 mL	135-100A



EARLY DETECTION • RAPID DIAGNOSIS • BETTER OUTCOMES



**For more information or to  
place an order, contact:**



**A PBM Group Company**

71 Veronica Avenue • Somerset • New Jersey • 08873

800.526.2125 • 732.246.3366 • Fax: 732.246.0570

[www.lifesignmed.com](http://www.lifesignmed.com)