

Status Mono

Immunoassay for the qualitative detection of infectious mononucleosis IgM heterophile antibodies in serum, plasma or whole blood



- **CLIA WAIVED** for whole blood
- Superior 99.9% sensitivity for early IM detection, specificity of 98.8%
- FDA cleared with no age restrictions
- Single reagent procedure provides results in 8 minutes
- External QC positive and negative controls included (for Items #68364 & 84M30)
- Convenient plastic sample collection pipettes

CPT CODE	
Whole Blood	86308QW
Plasma, Serum	86308

ORDERING INFORMATION		
PRODUCT	# OF TESTS	ITEM #
Status Mono (whole blood, serum or plasma)	25	68364
Status Mono (whole blood, serum or plasma)	30	84M30
Status Mono (whole blood only)	30	84W30
Status Mono (whole blood only)	10	84W10
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® Mono Controls	2 x 1 mL	200014

For more information or to place an order, contact: