



StatusTM

lifeSign

COVID-19/Flu A&B

A Rapid Immunoassay for the Simultaneous
Direct Detection and Differential Diagnosis of SARS-CoV-2,
Influenza Type A and Type B Antigen from nasopharyngeal swab specimens

Infections with these viruses may present similar symptoms.
Can you tell them apart? **We Can.**



Learn More

1 test + 1 specimen = 3 Results

The additive power of **Status**TM COVID-19/Flu A&B

- **Sensitivity and Specificity**
 - COVID-19 - Sensitivity 93.9%, Specificity 100%
 - Flu A - Sensitivity 91.4%, Specificity 95.7%
 - Flu B Sensitivity 87.6%, Specificity 95.9%
- **FDA Emergency Use Authorization (EUA)**
- **Visually read in 15 minutes**
- **Flocked nasopharyngeal swab for superior specimen collection and patient comfort**

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
Status TM Covid-19/Flu A&B	25	33225

Proudly made in the USA

For more information or to place an order, contact:

lifeSign

A PBM Group Company

85 Orchard Road • Skillman • New Jersey • 08558
800.526.2125 • 732.246.3366 • Fax: 732.246.0570 • www.lifesignmed.com