

Status DS

One Step Immunoassay For The Qualitative Detection Of
Drugs Of Abuse And Their Metabolites In Urine



- Accurate results with a 97-99% correlation to GC/MS
- Fast turnaround time with results in as little as 5 minutes
- Flexible combination of drug panels from a single drug up to 10 drug panels
- Test menu includes:

- | | |
|----------------------------------|------------------------|
| -Tricyclic Antidepressants (TCA) | -Opiates (OPI) |
| -Methadone (MTD) | -Phencyclidine (PCP) |
| -Amphetamine (AMP) | -Barbiturates (BAR) |
| -Methamphetamine (MET) | -Benzodiazepines (BZO) |
| -Marijuana (THC) | -Cocaine (COC) |